

MacroMedics®

OmniBoard™

The OmniBoard™ is the most versatile and advanced Modular Treatment Board of the new MacroMedics® Modular Patient Positioning Platform.


The OmniBoard offers unprecedented flexibility and accuracy in patient set-up. The OmniBoard Base can remain on the treatment couch during daily treatment and ensures precise patient positioning for all treatment areas. Combined with the modular MacroMedics positioning products, the OmniBoard becomes a dedicated device for head- and neck, lung, breast and thorax, liver and lower abdomen as well as SBRT treatments.

The unique detachable Breast Treatment Module offers maximum versatility. It can be easily attached to the OmniBoard Base within seconds and can be set in 3 treatment angles. The module is extremely light weight and has a comfortable 3D anatomical contoured design with optimized dosimetric properties. A detachable indexed bottom stop provides a rigid and secure patient support. The arm supports for breast, thorax and lung treatment are quickly transferable between the OmniBoard Base and Breast Treatment Module.

By adding the SBRT Secure Arches, the OmniBoard can be transformed to a state of the art SBRT Patient Positioning Platform for precise and accurate stereotactic treatment. The SBRT Secure Arches are available in 4 different sizes and can be customized with a respiratory suppression plate or vacuum cushions.

The OmniBoard is the first genuine “All in One” patient positioning solution available for radiotherapy treatment!

Contact MacroMedics or your local distributor for detailed information or an on-site product demonstration.



macromedics

SOLUTIONS IN RADIOTHERAPY

MacroMedics BV
Kouwe Hoek 18
2741 PX Waddinxveen
The Netherlands

Tel: +31(0)182 389777
Fax: +31(0)182 389778
info@macromedics.com
www.macromedics.com