





MR friendly

## Catheters

- Unprecedented Flexibility
- Kink Resistant
- Minimal Source friction
- Radio Opaque
- Sterile
- Sizes: 5F & 6F
- Treatment Length 300mm and 400mm
- Unique Flexi Stopper Buttons





- > 5 & 6 F 400mm Flexi Cath single Leader with button and one stopper button (sterile)
- > 5 & 6 F 400mm Flexi Cath single Leader (sterile)
- > 5 & 6 F 400mm Double Leader with two stopper buttons (sterile)
- > 5F & 6F 300mm & 400mm Flexi cath Single Leader with T-BALL and one stopper buttons (sterile)

**Catheters are suitable for all bodysites in brachytherapy and all remote afterloading systems. They are economic, safe, reliable and ideal for any clinic.**

New high-performance plastics have enabled the development of the best catheters in brachytherapy. They comply with the latest regulations on non-toxic medical implant materials. The low friction between source and catheter ensures high source dwell position accuracy

even in sharp curvatures. The catheters are easy to work with in daily practice. The precise, controlled diameter ensures a perfect fit if combined with buttons or templates.



Picture curtesy and T-Ball Flexi Cath designed by Dr. Rajendra Bhalavat



Kalyani Radiotherapy Specialty India (P) Ltd  
A-6 Vignesh Plaza,  
No. 11, First floor  
Thillainagar  
Trichy 620018  
Tamilnadu- India  
Tel: +91 431 2760433