

MacroMedics®

New: EagleBoard™

The new MacroMedics® EagleBoard™ offers quick and easy patient set-up for treatment of breast, lung and thorax.

A variety of indexed supports and grips for head and arms ensure comfortable as well as reproducible patient positioning for lung and thorax treatments. The tapered design for maximum CT bore clearance and the use of MR safe materials make the EagleBoard an ideal choice for all imaging modalities.


The New EagleBoard is a compact and cost-effective solution for quick and comfortable patient positioning.

The EagleBoard features:

- Tapered design for maximum CT bore clearance
- Easy adjustable positioning of the T- or U-Grip
- Grips can be fixed in two directions for a wide range of positions
- Full support of the upper arms reduces the need for additional devices such as vacuum cushions
- Adjustable and indexed position of MaxSupport 1 head support
- Compatible with the full range of MaxSupport head supports and blocks
- Short design allows a combination with additional thorax/abdominal immobilisation products such as respiratory suppression devices (EAMIS Lite)
- Indexable on all common couchtops for radiotherapy treatment and imaging
- Fully MR safe
- Comfortable and quick patient positioning



122800	EBTU	EagleBoard, 3P, MR safe including T- and U-Grip
122810	EBT	EagleBoard, 3P, MR safe including T-Grip
122820	EBU	EagleBoard, 3P, MR safe including U-Grip



macromedics

SOLUTIONS IN RADIOTHERAPY

MacroMedics BV
Kouwe Hoek 18
2741 PX Waddinxveen
The Netherlands

Tel: +31(0)182 389777
Fax: +31(0)182 389778
info@macromedics.com
www.macromedics.com