

MacroMedics® MultiBoard™

The MultiBoard™ is a versatile and accurate Modular Treatment Board and part of MacroMedics® new Modular Patient Positioning Platform.

This compact “All in One” patient positioning solution offers flexibility and allows quick, efficient and precise patient set-up for radiotherapy treatment.


The MultiBoard Base is anatomically shaped for patient comfort and can be indexed on the couchtop where it can remain during daily treatment. Its light weight carbon fiber structures are designed to offer optimized dosimetry.

Combined with the modular MacroMedics positioning products, the MultiBoard becomes a dedicated and accurate device for head and neck, lung and thorax, liver and lower abdomen as well as SBRT treatments.

Special thermoplastic masks for the different treatment areas can be attached to the MultiBoard. Additional fixation points in the shoulder area make it a dedicated solution for pediatric head and neck treatments. Other options are breast & thorax and pelvic masks.

By adding the SBRT Secure Arches, the MultiBoard can be transformed to a precise and comfortable high-end stereotactic device for accurate SBRT treatments. The Secure Arches are available in 4 different sizes and can be customized with a respiratory suppression plate or abdominal vacuum cushions.

Contact MacroMedics or your local distributor for detailed information or an on-site product demonstration.



macromedics

SOLUTIONS IN RADIOTHERAPY

MacroMedics BV
Kouwe Hoek 18
2741 PX Waddinxveen
The Netherlands

Tel: +31(0)182 389777
Fax: +31(0)182 389778
info@macromedics.com
www.macromedics.com